

A reliable and powerful source of compressed air

- ☑ Delivers dry, clean medical-grade compressed air
- ☑ Optimal cost-performance ratio
- ☑ Quiet performance suitable for use in any ICU

MEDICAL COMPRESSOR FOR VENTILATORS AND BUBBLE CPAPs

With its superior pressure and flow characteristics, the VENTIAIR compressor ensures a constant, clean, dry compressed air supply for medical ventilators. It can be used either as a primary air source or as an automatic backup to a hospital's air pipeline system for medical ventilators and NCPAP units.

TAURUS HEALTHCARE



VA11/ VA12



VA09



**VA09
FRONT**

Features

- Universally adaptable compressors
- Consistent performance in all conditions
- Simple and economical maintenance
- Low Power consumption
- Automatic water removal system



Compressor with digital display

Technical specifications	VA12	VA 11	VA 09
Nominal voltage and frequency	220 VAC+/- 10% 50/60Hz	220 VAC+/- 10% 50/60Hz	220 VAC+/- 10% 50/60Hz
Compressor efficiency	20 l.min-1/ 3,5 bar *	50 l.min-1/ 3,5 bar	50 l.min-1/ 3,5 bar
Air tank capacity	2 lit.	2 lit.	2 lit.
Filtration of air	5 µm	5 µm	5 µm
Dew point depression at 3bar, 20°C	5°C under the ambient temperature	5°C under the ambient temperature	5°C under the ambient temperature
Outlet connection	DISS 1160-A (3/4"-16 UNF)	DISS 1160-A (3/4"-16 UNF)	DISS 1160-A (3/4"-16 UNF)
Mode of the operation	Continuous	Continuous	Continuous
Automatic Back-up to central pipeline	on request	on request	on request
Noise level	≤49 dB [A]	≤49 dB [A]	≤49 dB [A]
Dimensions (l x b x h)	43.7 X 38.5 X 44 cm	43.7 X 38.5 X 44 cm	57 X 48 X 87 cm
Dimensions (l x b x h)-with castors	43.7 X 38.5 X 54.5 cm	43.7 X 38.5 X 54.5 cm	57 X 48 X 98 cm

* Bubble CPAP compressor



With Ventilator stand



Additional alarms
High Temp, low pressure



Circuit Support Arm

Compressors compatible with all ventilators and N-CPAP units



Taurus Healthcare Pvt. Ltd.

A5 A6 Pragati Chambers Ranjit Nagar
Commercial Complex, New Delhi-110008
Ph. : 011-45087820, Mobile : 9711203605, 9899966945
E-mail : aps087@gmail.com ; jaspreet@taurushealthcare.co.in
Reach us @: www.taurushealthcare.co.in

An Instromedix Group Company

