

Model - VA 11



A reliable and powerful source of compressed air

- Delivers dry, clean medical compressed air
- Optimal cost-performance ratio
- Quiet performance suitable for use in any ICU

MEDICAL COMPRESSOR FOR VENTILATORS AND BUBBLE CPAPS

With its superior pressure and flow characteristics the VENTIAIR compressor ensures a constant, clean, dry compressed air supply for medical ventilators. It can be used either as a primary air source or as an automatic backup to a hospital's air pipeline system for medical ventilators and NCPAP units.

FEATURES

- Consistent performance in all conditions
- Simple and economical maintenance
 - Low power consumption
- Automatic water removal system

Nominal voltage and frequency	220 VAC +/- 10% ; 50/60Hz
Compressor efficiency	50 L/min @ 3.5 Bar
Filtration of air	5 µm
Dew point depression at 40 L/min, 20°C	5°C under the ambient temperature
Outlet connection	Quick Coupler
Mode of the operation	Continuous
Automatic Back-up to central pipeline	On request
Noise level	≤50 dB
Dimensions (l x b x h)	44.5 x 43.8 x 44.0 cm
Dimensions (l x b x h)- with wheels	44.5 x 43.8 x 54.5 cm
Weight	46 kg

*Specifications may change without notice

TAURUS HEALTHCARE PVT. LTD.

(An ISO 13485 certified company)

B-26, Sector-57, Noida, UP-201301 • Mobile : 9711203605, 9899966945

E-mail : amarpal@taurushealthcareindia.com • Website : www.taurushealthcareindia.com

