



TAURUS VA 11 MEDICAL AIR COMPRESSOR

With its superior pressure and flow characteristics, the compressor ensures a constant, clean, dry compressed air supply.

It can be used either as a primary air source or as an automatic backup to a hospital's air pipeline system for respiratory devices.

FEATURES

- Consistent performance in all conditions
- Simple and economical maintenance
- Low power consumption
- Automatic water removal system
- Low noise

Nominal voltage and frequency	220 V~± 10%, 50 Hz ± 2% Hz
Output flow	45-50 LPM @ 50 PSI
Filtration of air	5 micron
Dew point depression at 3 bar, 20°C	3°C below the ambient temperature
Outlet connection	Quick Coupler
Mode of the operation	Continuous
Automatic Back-up to central pipeline	On request
Noise level	≤53 dB (A)
Dimensions (LBH)	45 x 44 x 44 cm
Dimensions (LBH)- with Castors	45 x 44 x 55 cm
Weight (with Castors)	47 kg

*Specifications may change without prior notice

DEVICE DESCRIPTION: Medical Air Compressor is used to supply clean dry oil free air to the respiratory devices. The device supports respiratory devices for air supply. The capacity is approximately 45-50 LPM at pressure range of 50 PSI. Compressor is intended to be operated by healthcare providers, physicians, nurses, and trained technicians.

USER ENVIRONMENT: Hospital (ICU/PICU/NICU) or sub-acute institution.

TAURUS HEALTHCARE PVT. LTD.

(An ISO 13485 certified company)

B-5, Sector-59, Noida, Gautam Buddha Nagar, Uttar Pradesh-201301, India
Mobile : 9711203605, 9899966945 • Phone : 011-250707558
E-mail : aps087@gmail.com, jaspreet@taurushealthcareindia.com • Website : www.taurushealthcareindia.com
Document No. THC/V1/PM-01 Rev 3.0 Dated 01.04.2025