

OXYGEN CONCENTRATOR



Oxygen Concentrator is intended to be used by the patients with respiratory disorders who require oxygen from external source. This device is used to provide oxygen only for medical use in home care or healthcare facility for providing oxygen therapy.

FEATURES

- Adjustable Flow Rate
- Single and dual flow meter option available
- Audio - Visual Alarms for Power Failure & Low - High Pressure (Optional)
- Provision for Nebulization & Oxygen Purity Indication (Optional)
- Light, and with Integrated Handle, Easy to Move, User Friendly.



SPECIFICATION	Model: Neumolite	Model: Neumolite 10
OXYGEN FLOW RANGE	0-5 L/min	0-10 L/min
OXYGEN CONCENTRATION	93% ±5%	93% ±5%
FLOWMETERS / OUTLET	1 or 2	1 or 2
MAXIMUM OUTLET PRESSURE	5-7 PSI	6-9 PSI
NOISE LEVEL	<50dB	<50dB
OPERATING TEMPERATURE RANGE	5°-32°C	5°-32°C
OPERATING HUMIDITY RANGE	20%-85% non-condensing	20%-85% non-condensing
POWER CONSUMPTION	350 W	500 W
VOLTAGE	220±10% V~; 50Hz	220±10% V~; 50Hz
WEIGHT	25 KG	29 KG

*Specifications may change without prior notice

TAURUS HEALTHCARE PVT. LTD.

(An ISO 13485 certified company)

B-5, Sector-59, Noida, Gautam Buddha Nagar, Uttar Pradesh-201301, India

Mobile : 9711203605, 9899966945 • Phone : 011-250707558

E-mail : aps087@gmail.com, jaspreet@taurushealthcareindia.com • Website : www.taurushealthcareindia.com

Document No. THC/V4/PM-01 Rev 3.0 Dated 01.04.2025