

TAURUS
T3



TAURUS
T5



Respiratory Heated Humidifier



- The Taurus T3/ Taurus T5 Respiratory Heated Humidifier is designed to warm and add humidity to gases delivered to patients requiring mechanical ventilation or other positive pressure breathing systems.
- Easy to use basic humidifier.
- The heater plate is controlled to a constant temperature which is determined by the setting of the smooth and easy to use heater control knob.
- Heater plate over temperature cut out for added patient safety.

Technical Specifications

Dimensions	135 mm X 170 mm X 155 mm
Weight	3 kg (without Chamber)
Display	3 digit, 14 mm LED Display (7 Segment) (Taurus T5 only)
Supply Voltage	230 V~±10%
Supply Frequency	50/60 Hz
Heater Control	1 - 9 (Corresponds to heater plate temperature of 45 - 75°C)
Heater Plate	90 W
Operating Temperature Range	18°C~26°C
Operating Humidity Range	< 85%
Heater Plate Over Temperature Cut Out	90±3°C
Digital Display of Patient Port Temperature	Taurus T5 Model only
High/Low Temperature Alarm	Taurus T5 Model only

TAURUS HEALTHCARE PVT. LTD.

(An ISO 13485 certified company)

B-5, Sector-59, Noida, Gautam Buddha Nagar, Uttar Pradesh-201301, India

Mobile : 9711203605, 9899966945 • Phone : 011-250707558

E-mail : aps087@gmail.com, jaspreet@taurushealthcareindia.com • Website : www.taurushealthcareindia.com

Document No. THC/V5/PM-02 Rev 3.0 Dated 01.04.2025